

ClearPeak One Step qRT-PCR Probe Kit

100 rxn

Cat No: 1SQP100

Shipping : Ship with blue ice.

Storage : Store at -20°C, protected from light.

For frequent use, the product may be stored at 2–8 °C. Avoid repeated freeze–thaw cycles.

General Information

ClearPeak One Step qRT-PCR Probe Kit is a one-step Real-Time RT-qPCR kit designed for rapid and sensitive detection of RNA targets.

Reverse transcription and quantitative PCR are performed consecutively within a single reaction mixture, eliminating the need for tube opening or additional reagent handling. This closed-tube format minimizes contamination risk and enhances workflow efficiency.

ClearPeak One Step qRT-PCR Probe Kit includes an advanced reverse transcriptase lacking RNase H activity, thereby preserving RNA integrity during cDNA synthesis. The enzyme exhibits high efficiency even with low RNA input and demonstrates strong performance with templates containing high GC content or complex secondary structures.

To improve specificity, the reverse transcriptase remains inactive at lower temperatures and is activated during the initial stages of thermal cycling. This feature helps reduce non-specific amplification caused by primer–dimer formation or non-specific primer binding, leading to improved accuracy in quantitative analysis.

The optimized buffer system enables both reverse transcription and PCR amplification to occur efficiently within the same reaction, maximizing overall performance. The kit offers high sensitivity, excellent specificity, and a broad dynamic range for precise quantification of target genes.

ROX reference dye is included to normalize fluorescence signal variations between wells. As different Real-Time PCR instruments require different ROX concentrations, users should select the appropriate setting based on their instrument specifications.

No ROX correction required: Roche LightCycler® 480, LightCycler® 96, Bio-Rad iCycler iQ, iQ5, CFX96

Low ROX required: ABI Prism® 7500/7500 Fast, QuantStudio® 3, 5, 6 Flex, 7 Flex, ViiA™ 7, Stratagene Mx3000/Mx3005P, Rotor-Gene® 3000

High ROX required: ABI Prism® 7000, 7300, 7700, 7900, ABI StepOne™ / StepOnePlus™, Eppendorf systems.

Instructions for using ROX:

ClearPeak 100x ROX Reference Dye should be added directly to the *ClearPeak* 2x One Step qRT-PCR Probe Buffer according to the ROX requirement of your real-time PCR instrument.

For 1 mL of *ClearPeak* 2x One Step qRT-PCR Probe Buffer, add:

High ROX instruments: add 20 µL of *ClearPeak* 100X ROX Reference Dye

Low ROX instruments: add 2 µL of *ClearPeak* 100X ROX Reference Dye

No ROX instruments: do not add ROX

After adding ROX, mix gently but thoroughly before use. Avoid repeated freeze–thaw cycles of the ROX-premixed buffer.

Protocol

1. Thaw the RNA template, primers, *ClearPeak* 2x One Step qRT-PCR Probe Buffer (ROX added according to the above protocol), *ClearPeak* One Step qRT-PCR Probe Enzyme Mix, and RNase-free water, and keep all components on ice until use.

Important Notes:

- Please gently mix the tube upside down before use and avoid foaming as much as possible. Avoid repeated freeze-thaw, as repeated freeze-thaw may cause a decrease in product performance.

2. Prepare the following reaction mixture

Component	Amount
<i>ClearPeak</i> 2x One Step qRT-PCR Probe Buffer	10 µl
<i>ClearPeak</i> One Step qRT-PCR Probe Enzyme Mix	1 µl
TaqMan Assay Probe (20x)	1 µl
RNA Template	<100ng
ddH ₂ O	up to 20 µl
Total	20 µl

3. Use the following reaction conditions to run the experiment. It is recommended to use a two-step PCR reaction program. If good experimental results cannot be obtained due to the use of primers with lower T_m values, a three-step PCR amplification can be attempted.

Two Step

Step	Temperature	Time	Cycles
Reverse Transcription	50°C	15-30 min	
Initial Denaturation	94°C	2 min	
Denaturation	94°C	15 sec	35-45
Annealing/Extension	60°C	30 sec	

Three Step

Step	Temperature	Time	Cycles
Reverse Transcription	50°C	15-30 min	
Initial Denaturation	94°C	2 min	
Denaturation	94°C	15 sec	35-45
Annealing	60°C	15 sec	
Extension	72°C	30 sec	

Important Notes:

- In practical operation, corresponding improvements and optimizations should be made to PCR Setup and PCR Reaction Conditions based on different templates, primer structures, and target fragment sizes.