

# ***ClearBand* 20x Tris-Buffered Saline, pH:7.6**

500 ml, 1000 ml

**Cat No:** TB2005, TB2010

**Shipping** : Ship at ambient temperature.

**Storage** : Store at room temperature. *ClearBand* 20x Tris-Buffered Saline is stable for 24 months. Use within 1 week once it is diluted to 1x.

## **General Information**

*ClearBand* 20x Tris-Buffered Saline (TBS) is a space-saving balanced pH stabilizing salt solution used for especially western blotting and ELISA procedures. *ClearBand* TBS enables washing without disturbing antibody-antigen binding interactions.

*ClearBand* 20x TBS should be diluted to 1x working solution with a pH  $7.6\pm 0.1$ .

*ClearBand* 20x TBS is prepared with ultra pure water and filter sterilized.