

NutriCulture DMEM, Low Glucose

500 ml

Cat No: DMEM-LG

Shipping : Ship at ambient temperature.

Storage : Store at 4-8°C with protection from light. Once opened, use within 6-8 weeks.

Formulation

Components	Concentration mg/L
Amino Acids:	
Glycine	30
L-Arginine HCl	84
L-Cystine 2HCl	62.57
L-Glutamine	584
L-Histidine HCl H ₂ O	42
L-Isoleucine	105
L-Leucine	105
L-Lysine HCl	146
L-Methionine	30
L-Phenylalanine	66
L-Serine	42
L-Threonine	95
L-Tryptophan	16
L-Tyrosine 2 Na	103.79
L-Valine	94

Components	Concentration mg/L
Vitamins:	
Choline Chloride	4
D-Calcium Pantothenate	4
Folic Acid	4
myo-Inositol	7.2
Nicotinamide	4
Pyridoxal HCl	4
Riboflavin	0.4
Thiamine HCl	4
Inorganic Salts:	
CaCl ₂ 2 H ₂ O	265
Fe(NO ₃) ₃ 9 H ₂ O	0.1
KCl	400
MgSO ₄	97.72
NaCl	6400
NaHCO ₃	3700
Na ₂ HPO ₄	109
Other Components:	
D-Glucose	1000
Sodium Pyruvate	110
Phenol Red Sodium Salt	15.9

Precautions and Disclaimer

This product is for research use only.