

***ClearBand* Acrylamide/Bisacrylamide**

37.5:1 %30 Ready-to-use Solution

100 ml, 500 ml

Cat No: AB3751-100, AB3751-500

Shipping : Ship at ambient temperature.

Storage : Store between 2°C and 8°C.

General Information

ClearBand Acrylamide/Bisacrylamide, 37.5:1 %30 Ready-to-use Solution is prepared using molecular biology grade acrylamide and bisacrylamide in ultrapure water. *ClearBand* Acrylamide/Bisacrylamide, 37.5:1 %30 Ready-to-use Solution is 0.22 µm filter sterilized and suitable for electrophoresis of both proteins and nucleic acids. It can be used for preparation of stacking and resolving gels during sodium dodecyl sulfate polyacrylamide gel electrophoresis.