

***ClearBand* 10x Phosphate Buffered Saline with Tween[®] 20, pH:7.4**

500 ml, 1000 ml

Cat No: PBST1005, PBST1010

Shipping : Ship at ambient temperature.

Storage : Store at room temperature. Use within 1 week once it is diluted to 1x.

General Information

ClearBand 10x Phosphate-Buffered Saline with Tween[®] 20 (PBST) is a balanced salt solution which is used for especially western blotting and ELISA procedures. *ClearBand* PBST enables washing without disturbing antibody-antigen binding interactions. *ClearBand* PBST does not contain calcium and magnesium.

ClearBand 10x PBST should be diluted to 1x working solution. 1x formulation contains 0.05% Tween 20.

ClearBand 10x PBST is prepared with ultra pure water and filter sterilized.