

# **ClearBand** Acrylamide/Bisacrylamide

## **37.5:1 %30 Ready-to-use Solution**

100 ml, 5x100 ml

**Cat No:** AB3751-100, AB3751-1005

**Shipping** : Ship at ambient temperature.

**Storage** : Store between 2°C and 8°C.

## **General Information**

**ClearBand** Acrylamide/Bisacrylamide, 37.5:1 %30 Ready-to-use Solution is prepared using molecular biology grade acrylamide and bisacrylamide in ultrapure water. **ClearBand** Acrylamide/Bisacrylamide, 37.5:1 %30 Ready-to-use Solution is 0.22 µm filter sterilized and suitable for electrophoresis of both proteins and nucleic acids. It can be used for preparation of stacking and resolving gels during sodium dodecyl sulfate polyacrylamide gel electrophoresis.