

NutriCulture Antibiotic Antimycotic Solution (100×)

100 ml

Cat No: AA100

Shipping : Ship with dry ice.

Storage : Store at -20°C.

General Information

NutriCulture Antibiotic Antimycotic Solution (100x) is an antibiotic mixture that can be effectively used to prevent gram-positive and gram-negative bacterial and fungal contamination in cell culture media. NutriCulture Antibiotic Antimycotic Solution is provided as a ready to use 100x concentrate and recommended final concentration is 10 ml/L.

Components

Penicillin G Sodium (10^7 Units/L); Streptomycin Sulfate (10000 mg/L); Amphotericin B (25 µg/ml); NaCl (9000 mg/L).