

EcoZyme RNase A

50 mg, 100 mg, 500 mg

Cat No: EZRA50, EZRA100, EZRA500

Shipping : Ship at ambient temperature.

Storage : Store at -20°C.

General Information

EcoZyme RNase A, from bovine pancreas, is a stable and highly reactive, DNase and protease-free, endoribonuclease that specifically degrades single stranded RNA at C and U residues. **EcoZyme** RNase A is pretreated to remove DNase I and is suitable for RNA digestion in plasmid and genomic DNA purification procedures.