

***ClearBand* 20x Tris-Buffered Saline with Tween® 20, pH:7.6**

500 ml, 1000 ml

Cat No: TBST2005, TBST2010

Shipping : Ship at ambient temperature.

Storage : Store at room temperature. Use within 1 week once it is diluted to 1x.

General Information

ClearBand 20x Tris-Buffered Saline with Tween® 20 (TBST) is a space-saving balanced pH stabilizing salt solution used for especially western blotting and ELISA procedures. ***ClearBand*** TBS enables washing without disturbing antibody-antigen binding interactions.

ClearBand 20x TBST should be diluted to 1x working solution. 1x formulation contains 0.1% Tween 20.

ClearBand 20x TBST is prepared with ultra pure water and filter sterilized.