

***ClearBand* Nuclease Free Ultra-Pure Water**

30 ml, 100 ml, 500 ml, 1000 ml, 10x30 ml, 10x100 ml, 5x500 ml

Cat No: DW003L, DW010L, DW05L, DW1L, DW003L10, DW010L10, DW05L5

Shipping : Ship at ambient temperature.

Storage : Store at room temperature.

General Information

ClearBand Nuclease Free Ultra-Pure Water is prepared under stringent conditions and is suitable for all applications in a molecular biology laboratory. *ClearBand* Nuclease Free Ultra-Pure Water goes through, deionization, reverse osmosis, UV-treatment, 0.2µm filtration, and double autoclave to ensure sterility, being free of nucleases and to completely eliminate the byproducts of DEPC.