

***NutriCulture* DMSO (Dimethyl Sulfoxide), Sterile**

Cell Culture Reagent

50 ml, 100 ml, 5x50 ml, 5x100 ml

Cat No: DMSO50, DMSO100, DMSO505, DMSO1005

Shipping : Ship at ambient temperature.

Storage : Store at room temperature. This product is stable for 24 months.

General Information

ClearBand DMSO (Dimethyl Sulfoxide), Sterile, is a colorless, polar, aprotic organic solvent used in chemical research. It is miscible with water and many other organic solvents.

ClearBand DMSO (Dimethyl Sulfoxide), Sterile, is also widely used in polymerase chain reactions (PCR) as a co-solvent, helping to inhibit the formation of secondary structures from DNA fragments and as a cryoprotectant vitrification agent for the preservation of cells, tissues and organs via prevention of ice crystal formation during cryopreservation.