

EcoZyme Proteinase K, powder

100 mg, 500 mg

Cat No: EZPK100, EZPK500

Shipping : Ship at ambient temperature.

Storage : Store at -20°C.

General Information

EcoZyme Proteinase K is a stable and highly reactive nonspecific serine protease that is useful for general digestion of proteins. **EcoZyme** Proteinase K is capable of inactivating RNases and DNases and is used in the isolation or preparation of high molecular weight nucleic acids.

EcoZyme Proteinase K is stable at various temperatures and may be used at 56°C for up to 4 hours, or 37°C for overnight incubations.