

***ClearBand* 1M Tris, pH 8.5**

500 ml, 1000 ml, 5x500 ml, 5x1000 ml

Cat No: T8505, T8510, T85055, T85105

Shipping : Ship at ambient temperature.

Storage : Store at room temperature.

General Information

ClearBand 1M Tris, pH 8.5 is a pre-mixed and pH-adjusted ready to use sterile-filtered solution. *ClearBand* 1M Tris, pH 8.5 can be diluted to desired concentration and can be used in molecular biology or general biochemistry applications.