

NutriCulture Trypsin-EDTA (0.25%)

100 ml

Cat No: TE025

Shipping : Ship with dry ice.

Storage : Store at -20°C.

General Information

NutriCulture Trypsin-EDTA (0.25%) contains EDTA and phenol red. NutriCulture Trypsin-EDTA (0.25%) is made from trypsin powder, which is a mixture of proteases derived from porcine pancreas. Trypsin-EDTA is widely used for cell culture applications including cell dissociation, routine cell culture passaging, and primary tissue dissociation.