

***NutriCulture* Trypsin-EDTA (0.05%)**

100 ml

Cat No: TE005

Shipping : Ship with dry ice.

Storage : Store at -20°C.

General Information

NutriCulture Trypsin-EDTA (0.05%) contains EDTA and phenol red. **NutriCulture** Trypsin-EDTA (0.05%) is made from trypsin powder, which is a mixture of proteases derived from porcine pancreas. Trypsin-EDTA is widely used for cell culture applications including cell dissociation, routine cell culture passaging, and primary tissue dissociation.