

***ClearBand* 20% SDS**

100 ml, 500 ml

Cat No: S100/20, S500/20

Shipping : Ship at ambient temperature.

Storage : Store at room temperature.

General Information

ClearBand 20% SDS (20% w/v) is sodium dodecyl sulfate prepared in 18 MΩ water. SDS is a detergent that is used for denaturation of proteins especially in polyacrylamide gel electrophoresis applications.

Precautions

If precipitation forms, incubate the solution at 50°C with occasional mixing until precipitate is completely dissolved.

Do not store ***ClearBand*** 20% SDS solution below 15°C.