

***NutriCulture* Fetal Bovine Serum**

50 ml

Cat No: FBS50

Shipping : Ship with dry ice.

Storage : Store below -20°C. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended.

General Information

NutriCulture FBS is especially rich in growth factors and is particularly low in antibodies, which may influence the cell culture work.

Quality control parameters of *NutriCulture* FBS include pH, osmolality, content of protein, albumin, IgG and hemoglobin, endotoxin level, sterility, mycoplasma detection and virus testing.

Country of Origin

South America

Precautions and Disclaimer

This product is for research use only.