

***ClearBand* Potassium Acetate**

3M, pH: 5.5

Ready-to-use Solution

100 ml

Cat No: PA100

Shipping : Ship at ambient temperature.

Storage : Store between 2-8°C.

General Information

ClearBand Potassium Acetate, 3M pH:5.5 Ready-to-use Solution is prepared using molecular biology grade potassium acetate in ultrapure water. *ClearBand* Potassium Acetate, 3M pH:5.5 Ready-to-use Solution is 0.22 µm filter sterilized and suitable for purification and precipitation of nucleic acids. It can be used for also precipitating dodecyl sulfate (DS) and DS-bound proteins allowing the removal of proteins from DNA.