

***ChemiPURE* Agarose, Electrophoresis Grade (low electroendosmosis)**

100 gr

Cat No: AGR01

Shipping : Ship at room temperature.

Storage : Store at room temperature (15–25°C) in a tightly closed container.

General Information

ChemiPURE Agarose, Electrophoresis Grade (low electroendosmosis), is suitable for routine DNA electrophoretic analysis. With good gel strength and clarity, it provides high purity for preparative and analytical electrophoresis of nucleic acids.

ChemiPURE Agarose is designed for laboratory nucleic acid electrophoresis workflows where clear gel formation and reliable separation are required.

Applications

- Routine DNA electrophoretic analysis
- Preparative electrophoresis of nucleic acids
- Analytical electrophoresis of nucleic acids

This product is sold for laboratory research use only.

CAS Number: 9012-36-6

Appearance: White powder

Storage Temperature: Room Temperature