

ChemiPURE Bovine Serum Albumin

5 gr, 10 gr, 25 gr, 50 gr

Cat No: BSA05, BSA10, BSA25, BSA50

Shipping : Ship with blue ice.
Storage : Store between 4-8°C.

General Information

ChemiPURE Bovine Serum Albumin (BSA) has a wide range of laboratory applications and is often used as a blocking agent in immune experiments, including ELISA and Western Blot. As carrier proteins, cross-linking them to haptens and other weak antigens can make them more immunogenic in antibody production. In restriction enzyme digestion reactions, BSA is often used as a reaction stabilizer for some enzymes and prevents it from adhering to tube walls and gun heads. BSA is also commonly used in the production process of biopharmaceuticals or as a nutrient in cell and microbial cultures. In addition, it is also used as a standard for protein quantitative detection. **ChemiPURE** Bovine Serum Albumin is prepared using high quality bovine plasma using Heat Shock method.

This product is sold for laboratory research use only.

CAS Number: 9048-46-8

Color: White to light yellow freeze-dried powder