

ChemiPURE Guanidine Hydrochloride

100 gr, 500 gr

Cat No: GH01, GH05

Shipping : Ship at ambient temperature.

Storage : Store at room temperature.

General Information

ChemiPURE Guanidine Hydrochloride is a strong chaotropic agent and a strong denaturant of proteins. It is an inhibitor of RNase. It can solubilize insoluble or denatured proteins such as inclusion bodies.

This product is sold for laboratory research use only.

CAS Number: 50-01-1

Molecular Weight: 95.5

Chemical Formula: $\text{CH}_5\text{N}_3\text{HCl}$

Purity: $\geq 99\%$

Color: White to off-white

Solubility: Water

Storage Temperature: Room Temperature