

***ClearBand* HRP Conjugated Goat anti-rabbit IgG Antibody**

500 µl

Cat No: HRPR5

Shipping : Ship with blue ice.

Storage : Store at -20°C with protection from prolonged exposure to light.

General Information

ClearBand HRP Conjugated Goat anti-rabbit IgG Antibody is labeled with horse radish peroxidase (HRP) and offer high sensitivity and performance for immunoblotting and ELISA applications.

ClearBand HRP Conjugated Goat anti-rabbit IgG Antibody reacts with whole molecule rabbit IgG based on immunoblotting and/or ELISA. It also reacts with the light chains of other rabbit immunoglobulins. This antibody has minimal cross-reaction with mouse, human and rat immunoglobulins, but it may cross-react with other immunoglobulins from other species.

Reactivity	Rabbit
Clonality	Polyclonal
Host	Goat
Applications	ELISA, WB
Conjugation	HRP
Formulation	Liquid in PBS containing 50% glycerol and 0.05% Proclin 300.
Purification	The antibody was purified by affinity chromatography.
Concentration	Lot-specific (please contact technical support for concentration and total µg amount)
Storage & Handling	The antibody solution should be stored undiluted at -20°C, and protected from prolonged exposure to light.