

***ClearBand* 20x Tris-Buffered Saline, pH:7.6**

500 ml, 1000 ml

Cat No: TB2005, TB2010

Shipping : Ship at ambient temperature.

Storage : Store at room temperature. *ClearBand* 20x Tris-Buffered Saline is stable for 24 months. Use within 1 week once it is diluted to 1x.

General Information

ClearBand 20x Tris-Buffered Saline (TBS) is a space-saving balanced pH stabilizing salt solution used for especially western blotting and ELISA procedures. *ClearBand* TBS enables washing without disturbing antibody-antigen binding interactions.

ClearBand 20x TBS should be diluted to 1x working solution with a pH 7.6 ± 0.1 .

ClearBand 20x TBS is prepared with ultra pure water and filter sterilized.