

***ClearBand* Nitrocellulose Membrane, 30 cm x 3 m, Roll**

0.22 µm, 0.45 µm

Cat No: NC022, NC045

Shipping : Ship at ambient temperature.

Storage : Store at room temperature.

General Information

ClearBand Nitrocellulose Membrane is dense 100% nitrocellulose combining the advantages of high protein binding capacity with low background and high membrane stability, which ensures easy handling and excellent signal-to-noise results.

ClearBand Nitrocellulose Membrane has excellent binding properties for western blotting, dot-blot assays, and other protein or nucleic acid methods. They are used for a wide range of molecular weight proteins and nucleic acids >300 bp.