

***ClearBand* dNTP Mix, 10 mM each**

1 ml, 5x1 ml

Cat No: ETP1, ETP5

Shipping : Ship with blue ice.

Storage : Store at -20°C for daily/weekly use. Store at -80°C for infrequent use (1-2 times per month). Avoid freeze and thaw cycles.

General Information

ClearBand dNTP (2'-deoxynucleoside 5'-triphosphate) Mix, 10 mM each consists of all four nucleotides (dATP, dCTP, dGTP, dTTP), each at a concentration of 10 mM, in water (pH: 7.4). *ClearBand* dNTP Mix is designed to save time and to provide higher reproducibility in PCR and other applications.