



ClearBand Fast Coomassie Blue Staining Solution

ready-to-use
solution



No toxicity: *ClearBand* Fast Coomassie Blue Staining Solution **does not contain toxic methanol and pungent odor of acetic acid** in conventional staining solution.

Quick: *ClearBand* Fast Coomassie Blue Staining Solution is formulated for a ready-to-use solution for quick staining (**as short as 1 minute**) of protein gels such as SDS-PAGE or non-denaturing PAGE.

To visit our website,
please scan the QR code

