


NutriCulture Mycoplasma Detection Kit



NutriCulture Mycoplasma Detection Kit offers a sensitive and specific option for the direct detection of mycoplasma in cell cultures via utilization of polymerase chain reaction (PCR).

The advantages of the *NutriCulture* Mycoplasma Detection Kit are its sensitivity, specificity, speed, cost efficiency, and the potential to screen a large number of samples.

Affordable: Low cost. An affordable alternative to competing kits.

! ⚠ *NutriCulture* Mycoplasma Detection Kit can be used together with *NutriCulture* Mycoplasma Elimination Kit, which is used for removal of mycoplasma from infected cell cultures with no cytotoxic side effects. *NutriCulture* Mycoplasma Elimination Kit can effectively eliminate mycoplasma strains, that are responsible for more than 85% of the contaminations in animal cell cultures, from experimentally contaminated and chronically infected cell lines.

To visit our website,
please scan the QR code

